

Don't leave it to chance be prepared with Coda[®].

The Coda[®] Inline[®] Filter patented technology is proven to effectively remove VOCs and CACs. They are CE certified and manufactured in a clean room. Each Coda[®] Inline[®] Filter contains a HEPA filter and can last between 3-6 months. All Coda[®] Inline[®] Filters are FDA 510(k) cleared, ISO 13485:2003 medical certified and made in the USA.

Coda[®] Xtra Inline[®] Filter
– Green with Aldasorb[®]



Pure activated carbon and Aldasorb[®]

Removes 98% of VOCs and aldehydes from gas

25 times absorptive capacity

Coda[®] Xtra Inline[®] Filter
– Blue



Pure activated carbon

20 times absorptive capacity

Coda[®] Regular Inline[®] Filter
– Purple



Pure activated carbon

10 times absorptive capacity

Smaller footprint

Reduced weight

