


Introducing our 4-Well GPS® Dish for IVF...

Our 4-well GPS dishes have been specifically designed and manufactured for IVF to meet the human ART needs for a high quality and easy to work 4-well dish.



Key Features and Big Benefits

- 4 large wells with GPS bottoms for easy viewing and manipulation.
- GPS bottom for embryo settling in the center of the wells for quick location of the embryo.
- Wells are designed to reduce shadows and have concise focal positioning.
- 4 large wells will hold up to 300 microliters of media and oil.
- 4 smaller outer wells hold up to 100 microliter of media and oil.
- Wide well design allows quick viewing and oil overlay of each well as well as the entire dish.
- No pyrogenics, clean medical grade.
- CE and ISO Registered, FDA 510(k) Cleared
- 1-cell MEA and < 0.03 LAL Tested

4-Well GPS®	4 large inner wells each with a 300 µl capacity, 4 smaller outer wells each with a 100 µl capacity. 10 dishes per sleeve	4GPS-010
µDrop GPS®	12 inner wells each with a 20 µl capacity, 32 mm inner ring holds 3 mL of oil, EmbryoAddress™ identification system, outer wells hold 4 mL (total) of oil or washing media. 10 dishes per sleeve.	MGPS-010
Universal GPS®	The 8 outer wells will each hold up to 100 µl of medium, and the 2 inner wells will each hold up to 150 µl. 10 dishes per sleeve.	UGPS-010
embryo GPS®	The 3 large center wells and 8 outer wells give you consistency. Large wells for flexibility in single or group culture. The 8 outer wells will each hold up to 50 µl of medium and the 3 inner wells will each hold up to 100 µl. Embryo culture has become simpler and safer. 10 dishes per sleeve.	EGPS-010
embryo corral®	The 2 large inner wells will each hold up to 120 µl (30 µl/quadrant) of medium and allow the flexibility in single or group culture. The 8 outer wells each will hold up to 50 µl of medium, and can be used for culturing, washing or preparation. 10 dishes per sleeve.	EMBC-010



IVFonline.com®

www.IVFOonline.com

T: 1-800-720-6375 ◆ F: 1-519-826-6947 ◆ Intl.: 001-519-826-5800 ◆ sales@LifeGlobal.com ◆ www.LifeGlobalGroup.com

©2015 LifeGlobal® Group, LLC

Dish_4-Well-2015-04-08-3478