

# LIFEGLOBAL®

## ...BEST QUALITY OILS FOR IVF CULTURE



### LifeGuard® Oil

- High viscosity
- For continuous culture to protect embryos from VOCs

100 ml LGUA-100  
500 ml LGUA-500  
100 ml AMLG-100



### LiteOil®

- Medium viscosity
- For culture & micromanipulation
- High Quality mineral oil

100 ml LGOL-100  
500 ml LGOL-500  
100 ml AMLO-100



### Paraffin Oil P.G.

- Low viscosity
- For culture & micromanipulation
- High Quality 100% paraffin oil

100 ml LGPO-100  
500 ml LGPO-500  
100 ml AMPO-100

## Protect your culture system with peroxide-free oil.

- High Quality & Consistency
- Sterile filtered
- Raw materials are washed with ultra-pure water to remove water soluble contaminants
- 1-cell MEA, Endotoxin Tested & cell count
- 24 Month Shelf Life
- CE Certified and ISO 13485 Certified
- Specifically designed & tested for IVF
- Pharmaceutical grade
- Maintain temperature & media osmolarity of the culture media
- Minimize pH shifts of the culture media
- < 0.001 Peroxide
- Store dark at 15-25°C
- Ready-to-use after equilibration at 37°C and 5-6% CO<sub>2</sub>