


LIFEGLOBAL®

...BEST QUALITY OILS FOR IVF CULTURE



LifeGuard® Oil

- High viscosity
 - For continuous culture to protect embryos from VOCs
- 100 ml LGUA-100
500 ml LGUA-500
100 ml AMLG-100



LiteOil®

- Medium viscosity
 - For culture & micromanipulation
 - High Quality mineral oil
- 100 ml LGOL-100
500 ml LGOL-500
100 ml AMLO-100



Paraffin Oil P.G.

- Low viscosity
 - For culture & micromanipulation
 - High Quality 100% paraffin oil
- 100 ml LGPO-100
500 ml LGPO-500
100 ml AMPO-100

Protect your culture system with peroxide-free oil.

- High Quality & Consistency
- Sterile filtered
- Raw materials are washed with ultra-pure water to remove water soluble contaminants
- 1-cell MEA, Endotoxin Tested & cell count
- 24 Month Shelf Life
- FDA Cleared & CE Marked
- Specifically designed & tested for IVF
- Pharmaceutical grade
- Maintain temperature & media osmolarity of the culture media
- Minimize pH shifts of the culture media
- < 0.001 Peroxide
- Store dark at 15-25°C
- Ready-to-use after equilibration at 37°C and 5-6% CO₂