

# GPS DISHWARE®

...A SPECIALLY DESIGNED DISH FOR EVERY IVF PROCEDURE



## Mini GPS® 38


- Ideal for embryo culture
  - 7 - 20 µl wells
- MI38-010



## 38 Special GPS®

- Ideal for embryo culture
  - 12 - 15 µl wells
  - 2 washing bays
- SP38-010

- GPS design for quick embryo location
- Enhanced optics and visibility
- Raised lids for rapid and constant equilibration
- CE Cleared and ISO Registered, FDA 510(k) Cleared
- 1-cell MEA, Endotoxin Tested (LAL), Sterile
- Non-toxic medical grade polystyrene
- Lower VOC & breathable packaging
- Ideal for PGD



## µDrop GPS®

- Ideal for continuous culture
  - 12 - 20 µl wells
- MGPS-010




## 4-Well GPS®

- Ideal for large volume culture & vitrification
  - 4 - 300 µl wells
  - 4 - 100 µl wells
- 4GPS-010



- Improved safety with GPS design



## embryo corral®

- Ideal for group culture
  - 2 - 120 µl wells
  - 8 - 50 µl wells
- EMBC-010



## embryo GPS®

- Ideal for PGD cases
  - 3 - 100 µl wells
  - 8 - 50 µl wells
- EGPS-010



## Universal GPS®

- Ideal for all procedures
  - 2 - 150 µl wells
  - 8 - 100 µl wells
- UGPS-010